***Echinococcus multilocularis* (Em18) IgG ELISA**

Imunoenzimski test za serološko praćenje humane alveolarne ehinokokoze

96 testova na pojedinačnim jamiicama za in vitro dijagnostičku upotrebu i profesionalnu laboratorijsku upotrebu

Uputstvo za upotrebu za proizvod N° 9310
UDI-DI: 07640158219317

Predviđena upotreba:

Bordier *Echinococcus multilocularis* (Em18) IgG ELISA kit je namenjen za kvalitativnu detekciju IgG antitela na antigen Em18 parazita *Echinococcus multilocularis* u humanom serumu. Ovaj test je namenjen za postoperativno i/ili postterapijsko praćenje inficiranih pacijenata.

Pozadina:

Alveolarna ehinokokoza (AE) je uzrokovana larvalnim stadijumom *Echinococcus multilocularis*, pantljičare koja se javlja kod lisica, kojota, pasa i drugih kanida. Ljudi se mogu zaraziti slučajnim unošenjem jaja pantljičare putem kontaminirane hrane ili vode. Larve *E. multilocularis* ne sazrevaju u potpunosti u plodne ciste kod ljudi, ali će kontinuirana proliferacija vezikula koje invadiraju i uništavaju ga nalik tumoru prouzrokovati disfunkciju jetre. Parazit se može proširiti na druge organe, kao što su pluća i mozak. Glavni simptomi su bol u stomaku, astenija, hepatomegalija i žutica. Dijagnoza se zasniva na tehnikama vizuelizacije kao što su CT skeniranje za vizuelnu detekciju parazitskih masa i odgovarajućih difuznih struktura nalik cistama. Serološki testovi se koriste za skrining rizične populacije i za praćenje AE pacijenata nakon lečenja.

Princip i prezentacija:

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Reference:

- Tappe, D., Frosch, M., Sako, Y., Itoh, S., Gruner, B., Reuter, S., et al. (2009) Close relationship between clinical regression and specific serology in the follow-up of patients with alveolar echinococcosis in different clinical stages. *Am. J. Trop. Med. Hyg.* **80** : 792-7.
- Ishikawa, Y., Sako, Y., Itoh, S., Ohtake, T., Kohgo, Y., Matsuno, T., et al. (2009) Serological monitoring of progression of alveolar echinococcosis with multiorgan involvement by use of recombinant Em18. *J. Clin. Microbiol.* **47** : 3191-6.
- Tappe, D., Sako, Y., Itoh, S., Frosch, M., Gruner, B., Kern, P., Ito, A. (2010) Immunoglobulin G subclass responses to recombinant Em18 in the follow-up of patients with alveolar echinococcosis in different clinical stage. *Clin Vaccine Immunol.* **17** : 944-8.
- Ammann, R.W., Stumpe, K.D.M., Grimm, F., Deplazes, P. et al. (2015) Outcome after discontinuing long-term benzimidazole treatment in 11 patients with non-resectable alveolar echinococcosis with negative FDG-PET/CT and anti-Em18/3-10 serology. *PLoS Negl Trop Dis.* **9**.
- Sulima, M., Szostakowska, B., Nahorsk, W., Sikorska, K., Wolyniec, W. and Waz, P. (2019) The usefulness of commercially available serological tests in the diagnosis and monitoring of treatment in patients with alveolar echinococcosis. *Clin. Exp. Hep.* **5** : 327-333.
- Hotz, J.F., Peter, L., Kapp-Schworer, S., Theis, F., Eberhardt, N., Essig, A. et al. (2022) Evaluation of Serological Markers in Alveolar Echinococcosis Emphasizing the Correlation of PET-CTI Tracer Uptake with RecEm18 and Echinococcus-Specific IgG. *Pathogens* **11**:20239.

BORDIER AFFINITY PRODUCTS S.A.

📍 Chemin de Chatanerie 2, 1023 Crissier, Switzerland.

☎ +41 21 633 31 67 ✉ cb@bordier.ch 🌐 www.bordier.ch